Material Safety Data Sheet

Date prepared: August, 29th, 2012 Date revised: N/A (first edition)

Section 1 - General Information

Chemical Name: Diphenyl phthalate

Synonyms: NA
CAS Reg. Number: 84-74-2
Catalog Number: 1ST1114
Manufacturer's Name: A Chemtek, Inc.

100 Barber Avenue

Worcester, MA 01606, USA

Telephone Number: 508-471-4121 **Fax Number:** 508-845-9201

For R & D use only

Section 2 - Hazardous Ingredients/Identity Information

Classification of the substance or mixture

GHS Classification in accordance with 29 CFR 1910 (OSHA HCS)

Reproductive toxicity (Category 1B), H360 Acute aquatic toxicity (Category 1), H400 Chronic aquatic toxicity (Category 1), H410

GHS Label elements, including precautionary statements



Pictogram
Signal word
Danger

Hazard statement(s)

H360 May damage fertility or the unborn child.

H410 Very toxic to aquatic life with long lasting effects.

Precautionary statement(s)

P201 Obtain special instructions before use.

P202 Do not handle until all safety precautions have been read and

understood.

P273 Avoid release to the environment.

P281 Use personal protective equipment as required.

P308 + P313 IF exposed or concerned: Get medical advice/ attention.

P391 Collect spillage. P405 Store locked up.

P501 Dispose of contents/ container to an approved waste disposal plant.

HMIS Classification

Health hazard: 1 Chronic Health Hazard: * Flammability: 1 Physical hazards: 0

NFPA Rating

Health hazard: 2 Fire: 1 Reactivity Hazard: 0

Section 3 - Composition/Information on Ingredients

 $\begin{array}{lll} Formula: & C_{16}H_{22}O_4 \\ Molecular Weight: & 278.34 \text{ g/mol} \\ CAS-No. & 84-74-2 \end{array}$

Section 4 - First Aid Measures

GENERAL ADVICE: Consult a physician. Show this safety data sheet to the doctor in attendance.

Move out of dangerous area.

Material Safety Data Sheet

Date prepared: August, 29th, 2012 Date revised: N/A (first edition)

EYE CONTACT: Rinse thoroughly with plenty of water for at least 15 minutes and consult a

physician.

SKIN CONTACT: Wash off with soap and plenty of water. Consult a physician.

INHALATION: If breathed in, move person into fresh air. If not breathing, give artificial

respiration. Consult a physician.

INGESTION: Never give anything by mouth to an unconscious person. Rinse mouth with

water. Consult a physician.

Section 5 - Fire-Fighting Measures

Flammable properties

Flash point not available Ignition temperature not available

Suitable extinguishing media

Use water spray, alcohol-resistant foam, dry chemical or carbon dioxide.

Special protective equipment for fire-fighters

Wear self-contained breathing apparatus for firefighting if necessary.

Section 6 - Accidental Release Measures

Personal precautions Use personal protective equipment. Avoid breathing vapours, mist

or gas. Ensure adequate ventilation. Evacuate personnel to safe areas.

Environmental precautions Prevent further leakage or spillage if safe to do so. Do not let

product enter drains.

Methods for cleaning up Sweep up and shovel. Keep in suitable, closed containers for

disposal.

Section 7 - Handling and Storage

Handling Avoid inhalation of vapour or mist.

Storage Keep container tightly closed in a dry and well-ventilated place. Containers which

are opened must be carefully resealed and kept upright to prevent leakage.

Store under inert gas. Moisture sensitive.

Section 8 - Exposure Controls/Personal Protection

Components with workplace control parameters

Personal protective equipment

Respiratory protection Where risk assessment shows air-purifying respirators are appropriate use a

full-face respirator with multi-purpose combination (US) or type ABEK (EN 14387) respirator cartridges as a backup to engineering controls. If the respirator is the sole means of protection, use a full-face supplied air respirator. Use respirators and components tested and approved under appropriate government standards such as NIOSH (US) or CEN (EU). Handle with gloves. Gloves must be inspected prior to use. Use proper

Handle with gloves. Gloves must be inspected prior to use. Use proper

glove removal technique (without touching glove's outer surface) to avoid skin contact with this product. Dispose of contaminated gloves after use in accordance with applicable laws and good laboratory practices. Wash and

dry hands.

Eye protection Safety glasses with side-shields conforming to EN166 Use equipment for

eye protection tested and approved under appropriate government standards

such as NIOSH (US) or EN 166(EU).

Skin protection Impervious clothing, The type of protective equipment must be selected

according to the concentration and amount of the dangerous substance at the

specific workplace.

Hygiene measures Handle in accordance with good industrial hygiene and safety practice.

Wash hands before breaks and at the end of workday.

Section 9 - Physical and Chemical Properties

Appearance: Form Liquid

Material Safety Data Sheet

Date prepared: August, 29th, 2012 Date revised: N/A (first edition)

Color colorless

Safety data

pH: not available Melting point: $-35.0 \,^{\circ}\text{C} \, (-31.0 \,^{\circ}\text{F})$ Boiling point: $340.0 \,^{\circ}\text{C} \, (644.0 \,^{\circ}\text{F})$

Flash point: $171.0 \,^{\circ}\text{C} \, (339.8 \,^{\circ}\text{F}) - \text{closed cup}$

Ignition temperature: not available

Lower explosion limit: Lower explosion limit: 0.47 %(V)

Upper explosion limit: not available Water solubility: not available

Section 10 - Stability and Reactivity

Storage stability Stable under recommended storage conditions.

Materials to avoid Strong oxidizing agents.

Hazardous decomposition products Hazardous decomposition products formed under fire conditions.

Carbon oxides.

Section 11 - Toxicological Information

Acute toxicity:

Oral LD50

LD50 Oral - rat - 8,000 mg/kg

Inhalation LC50

LC50 Inhalation - rat - 4,250 mg/m3

Dermal LD50

LD50 Dermal - rabbit - > 20,860 mg/kg Other information on acute toxicity

no data available

Irritation and corrosion: not available

Sensitization: not available **Respiratory or skin sensitization:** not available

Chronic exposure No component of this product present at levels greater than or equal to 0.1%

is identified as probable, possible or confirmed human carcinogen by IARC,

ACGIH, NTP, or OSHA.

Section 12 - Ecological Information

Elimination information (persistence and degradability)

not available

Ecotoxicity effects: not available

Further information on ecology: not available

Section 13 - Disposal Considerations

Product Offer surplus and non-recyclable solutions to a licensed disposal company.

Contact a licensed professional waste disposal service to dispose of this

material.

Contaminated packaging: Dispose of as unused product.

Section 14- Transport Information

DOT (US):

UN number: 3082 Class: 9 Packing group: III

Proper shipping name: Environmentally hazardous substances, liquid, n.o.s. (Benzyl butyl phthalate)

Reportable Quantity (RQ): 10 lbs

Marine pollutant: No

Poison Inhalation Hazard: No

IMDG:

UN number: 3082 Class: 9 Packing group: III EMS-No: F-A, S-F

Material Safety Data Sheet

Date prepared: August, 29th, 2012 Date revised: N/A (first edition)

Proper shipping name: ENVIRONMENTALLY HAZARDOUS SUBSTANCE, LIQUID, N.O.S.

(Benzyl butyl phthalate)

Marine pollutant: Marine pollutant

IATA:

UN number: 3082 Class: 9 Packing group: III

Proper shipping name: Environmentally hazardous substance, liquid, n.o.s. (Benzyl butyl phthalate)

Section 15 - Regulatory Information

SARA 302 Components

SARA 302: No chemicals in this material are subject to the reporting requirements of SARA Title III, Section 302.

SARA 313 Components

The following components are subject to reporting levels established by SARA Title III, Section 313:

Dibutyl phthalate CAS-No. 84-74-2 Revision Date 2007-07-01

SARA 311/312 Hazards

Chronic Health Hazard

Massachusetts Right To Know Components

No components are subject to the Massachusetts Right to Know Act.

Pennsylvania Right To Know Components

Dibutyl phthalate CAS-No. 84-74-2 Revision Date 2007-07-01

New Jersey Right To Know Components

Dibutyl phthalate CAS-No. 84-74-2 Revision Date 2007-07-01

California Prop. 65 Components

WARNING: This product contains a chemical known to the State of California to cause birth defects or

other reproductive harm.

Dibutyl phthalate CAS-No. 84-74-2 Revision Date 2008-07-01

Section 16 - Other Information

The above information is believed to be correct but does not purport to be all-inclusive and shall be used only as a guide. A Chemtek shall not be held liable for any damage resulting from handling or from contact with the above product.