Material Safety Data Sheet

Date revised: N/A (first edition)

Section 1 - General Information	
Chemical Name:	Sulfacetamide
CAS Reg. Number:	144-80-9
Catalog Number:	1ST4001
Manufacturer's Name:	A Chemtek, Inc.
	100 Barber Avenue
	Worcester, MA 01606, USA
Telephone Number:	508-471-4121
Fax Number:	508-845-9201
For R & D use only	

Section 2 - Hazardous Ingredients/Identity Information

Classification of the substance or mixture

Not a hazardous substance or mixture according to Regulation (EC) No. 1272/2008.

Not a hazardous substance or mixture according to EC-directives 67/548/EEC or 1999/45/EC.

Label elements

The product does not need to be labelled in accordance with EC directives or respective national laws. Other hazards - none

Section 3 - Composition/Information on Ingredients

Formula:	$C_8H_{10}N_2O_3S$
Molecular Weight:	214.24g/mol
CAS-No.	144-80-9

Section 4 - First Aid Measures

If inhaled If breathed in, move person into fresh air. If not breathing give artificial respiration. In case of skin contact Wash off with soap and plenty of water. In case of eye contact Flush eyes with water as a precaution. If swallowed Never give anything by mouth to an unconscious person. Rinse mouth with water.

Section 5 - Fire-Fighting Measures

Suitable extinguishing media Use water spray, alcohol-resistant foam, dry chemical or carbon dioxide. Special protective equipment for fire-fighters Carbon oxides, nitrogen oxides (NOx), Sulphur oxides.

Section 6 - Accidental Release Measures

Personal precautions	Avoid dust formation.
Environmental precautions	Do not let product enter drains.
Methods for cleaning up	Sweep up and shovel. Keep in suitable, closed containers for disposal.

Section 7 - Handling and Storage

Provide appropriate exhaust ventilation at places where dust is formed. Normal measures for Handling preventive fire protection. Keep container tightly closed in a dry and well-ventilated place. Storage

Section 8 - Exposure Controls/Personal Protection

Control parameters Components with workplace control parameters **Exposure controls Appropriate engineering controls** General industrial hygiene practice. Personal protective equipment

MSDS

Material Safety Data Sheet

Eye/face protection

Use equipment for eye protection tested and approved under appropriate government standards such as NIOSH (US) or EN 166(EU).

Skin protection

Handle with gloves. Gloves must be inspected prior to use. Use proper glove removal technique (without touching glove's outer surface) to avoid skin contact with this product. Dispose of contaminated gloves after use in accordance with applicable laws and good laboratory practices.

Wash and dry hands.

The selected protective gloves have to satisfy the specifications of EU Directive 89/686/EEC and the standard EN 374 derived from it.

Body Protection

Choose body protection in relation to its type, to the concentration and amount of dangerous substances, and to the specific work-place., The type of protective equipment must be selected according to the concentration and amount of the dangerous substance at the specific workplace.

Respiratory protection

Respiratory protection is not required. Where protection from nuisance levels of dusts are desired, use type N95 (US) or type P1 (EN 143) dust masks. Use respirators and components tested and approved under appropriate government standards such as NIOSH (US) or CEN (EU).

Section 9 - Physical and Chemical Properties

Appearance:	
Form	Powder
Color	beige
Safety data	
pH:	not available
Melting point:	Melting point/range: 182-184 °C
Boiling point:	not available
Flash point:	not available
Ignition temperature:	not available
Lower explosion limit:	not available
Upper explosion limit:	not available
Water solubility:	not available
Other safety information	Bulk density 0.5g/l

Section 10 - Stability and Reactivity

Storage stability	Stable under recommended storage conditions
Materials to avoid	Strong oxidizing agents
Hazardous decomposition products	Other decomposition products - no data available

Section 11 - Toxicological Information

Acute toxicity:	LD50 Oral - mouse - 16.500 mg/kg
Irritation and corrosion:	not available
Sensitization:	not available
Respiratory or skin sensitization:	Prolonged or repeated exposure may cause allergic reactions in certain
	sensitive individuals.

Chronic exposure

No component of this product present at levels greater than or equal to 0.1% is identified as probable, possible or confirmed human carcinogen by IARC.

Potential health effects

Inhalation:	May be harmful if inhaled. May cause respiratory tract irritation.
Skin:	May be harmful if absorbed through skin. May cause skin irritation.
Eyes:	May cause eye irritation.
Ingestion:	May be harmful if swallowed.

Section 12 - Ecological Information

Elimination information (persistence and degradability)

MSDS

Date prepared: May 28, 2012

Material Safety Data Sheet

Date revised: N/A (first edition)

	not available
Ecotoxicity effects:	not available
Further information on ecology:	not available

Section 13 - Disposal Considerations

Product Contaminated packaging Offer surplus and non-recyclable solutions to a licensed disposal company. Dispose of as unused product.

Section 14- Transport Information

DOT (US):	Not dangerous goods
IMDG:	Not dangerous goods
IATA:	Not dangerous goods

Section 15 - Regulatory Information

This safety datasheet complies with the requirements of Regulation (EC) No. 1907/2006.

Section 16 - Other Information

The above information is believed to be correct but does not purport to be all-inclusive and shall be used only as a guide. A Chemtek shall not be held liable for any damage resulting from handling or from contact with the above product.